


# Exemplo | Examples

- Exemplo 1 | Example #1

Paciente masculino, 96 anos, com síndrome coronariana aguda. Em uso de carvedilol e espironolactona. Ao ECG, apresenta ritmo de marca-passo artificial, sem alterações de sensibilidade ou captura, compatível com sistema de estimulação cardíaca normofuncionante.


*Male patient, 96 years old, presents with acute coronary syndrome. He takes carvedilol and spironolactone daily. The ECG shows a rhythm from an artificial pacemaker, with stable sensitivity and capture, indicating a normally functioning cardiac stimulation system.*



- Exemplo 2 | Example #2

Paciente feminina, 68 anos, com queixa de epigastralgia. Tem diagnóstico de hipertensão arterial sistêmica. Em uso de betabloqueador e diurético. Ao ECG, apresenta ritmo de marca-passo artificial, sem alterações de sensibilidade ou captura, compatível com sistema de estimulação cardíaca normofuncionante.

*A 68-year-old female patient is experiencing epigastric pain. She has been diagnosed with systemic arterial hypertension and takes beta blockers and diuretics on a daily basis. An ECG reveals an artificial pacemaker rhythm. There are no changes in sensitivity or capture, indicating that the cardiac stimulation system is functioning normally.*



---

Revision #1

Created 22 May 2025 16:24:40 by Gabriela Miana

Updated 22 May 2025 16:28:24 by Gabriela Miana